90-46450-01

July 2010

RESTYLANE-L™ INSTRUCTIONS FOR USE

Formula

Contents:
Hyaluronic acid, stabilized 20 mg/mL
Lidocaine hydrochloride 0.3%
Phys. sodium chloride solution, pH 7 q.s.

Description

Restylane-L is a clear, transparent and viscous gel supplied in a glass syringe together with a 29 G Thin Wall (TW) or 30 G needle(s). The product is for single use only. *Restylane-L* is a unique form of non-animal, stabilized hyaluronic acid (NASHA™). Hyaluronic acid is a natural polysaccharide which occurs as an important structural element in the skin and in subcutaneous and connective tissues as well as in the synovial tissue and fluid. Hyaluronic acid belongs to a group of very few substances which are identical in all living organisms.

Mode of action

Restylane-L acts by adding volume to the tissue, thereby restoring the skin contours or enhancing the lips to the desired level of correction. *Restylane-L* is naturally integrated into the tissue and will in time undergo isovolemic degradation.

Restylane-L has been enhanced with the addition of Lidocaine in order to reduce patient discomfort during treatment.

Indication and usage

Restylane-L is intended to be used for facial tissue augmentation. It is recommended that the product be used for the correction of wrinkles and for lip enhancement. It should be injected into the middle part of the dermis layer of the facial skin. With cutaneous contour deformities the best results are obtained if the defect can be manually stretched to the point where it is eliminated. The degree and duration of the correction depend on the character of the defect treated, the tissue stress at the implant site, the depth of the implant in the tissue and the injection technique. Markedly indented defects may be difficult to correct. For the correction of thin superficial lines *Restylane Fine Lines*® is recommended. For lip enhancement both *Restylane*® and *Perlane*® can be used. For shaping the contours of the face and for the correction of folds *Perlane* is recommended. Also combinations of *Restylane*, *Perlane* and *Restylane Fine Lines* can be used. Please consult the *Perlane* and the *Restylane Fine Lines* instructions for Use for more information.

Contraindications

Do not use *Restylane-L* on patients with known hypersensitivity to lidocaine or amide type local anaesthetics.

Warning

Restylane-L is only intended for use as an intradermal implant. Do not reutilize *Restylane-L*. Do not inject intravascularly. Do not mix with other products.

Precautions

Normal precautions associated with intradermal injections must be observed. Like any such procedure, the implantation of *Restylane-L* is associated with an inherent risk of infection. *Restylane-L* should not be used in or near anatomic sites where there is active skin disease, inflammation or related conditions. Do not use *Restylane-L* together with any other injectable implant, except for *Restylane Fine Lines* and *Perlane* products. *Restylane-L* should not be injected into an area where a permanent implant has been placed.

Restylane-L should not be used for patients with unattainable expectations. The patient should be informed that he or she should not expose the treated area to intense heat (e.g. solarium and sunbathing) or extreme cold at least until the initial swelling and redness have resolved.

If laser treatment, chemical peeling or any other procedure based on active dermal response is considered after treatment with *Restylane-L*, there is a theoretical risk of eliciting an inflammatory reaction at the implant site. This also applies if *Restylane-L* is administered before the skin has healed completely after such a procedure.

There is a potential risk with the procedure that the material could be inadvertently injected into dermal blood vessels, which could lead to vascular occlusion in an end-artery with corresponding consequences. *Restylane-L* has not been tested in pregnant or lactating women or in children.

Anticipated side-effects

After the injection of *Restylane-L*, some common injection-related reactions might occur. These reactions include erythema, swelling, pain, itching, discoloration or tenderness at the implant site. Typically resolution is spontaneous within a few days after injection into the skin and within a week after injection into the lips. Additionally, temporary palpable lumpiness has been noted after the use of *Restylane-L* in some patients. The composition of *Restylane-L* is set to balance a normal tissue pressure. However, because the tissue pressure is sometimes disturbed to a higher value such as during oedema or to a lower value such as during dehydration, a small but significant change (swelling or shrinkage) may occur.

Adverse events

Reactions thought to be of a hypersensitivity nature have been reported in about 1 in every 10,000 treatments. These have consisted of swelling and induration at the implant site, sometimes with oedema in the surrounding tissues. Erythema, tenderness and rarely acneiform papules may also occur. The reactions have started either shortly after injection or after a delay of 2-4 weeks and have generally been described as mild to moderate and self-limiting, with an average duration of 2 weeks. In pronounced cases a short course of oral corticosteroids may prove effective. In addition, rare cases (less than 1 in 15,000 treatments) of granuloma formation, superficial necrosis and urticaria have been reported. Patients who have experienced this type of reaction should not be retreated with a *Restylane* product.

Adverse events must be reported to Medicis Aesthetics Canada Ltd.

Interactions

Treatment with *Restylane-L* in combination with other drugs and devices has not been tested.

Assembly of 30 G needle to syringe

For safe use of *Restylane-L* it is important that the needle is properly assembled.

See pictures A through D.

- Unscrew the tip cap (A) of the syringe carefully.
- Take a loose grip on the narrow part of the needle shield and mount the needle on the Luer-lok® (B) by screwing until you feel some counterpressure.
- Take a new firm grip on the wider part of the needle shield (C). **Press and turn it a further 90°** (a quarter of a turn).
- Pull off the needle shield.

A.

B.

C.


D.

Assembly of 29 G (TW) needle to syringe

Use the thumb and forefinger to hold firmly around both the glass syringe barrel and the Luer-Lok adapter. Grasp the needle shield with the other hand. To facilitate proper assembly, both push and rotate firmly.

Dosage and administration

Before the treatment, the patient's suitability for the treatment and the need for pain relief should be assessed. Normally, no anaesthesia is necessary when treating wrinkles. For lip augmentation, anaesthesia through a nerve block can be used. The patient should be informed about the indications, expected result, contraindications, precautions, warnings and potential adverse events. The treatment site should be cleaned with a suitable antiseptic solution. *Restylane-L* is administered using a 29 G thin wall or 30 G needle by injecting the material into the dermis. If *Restylane-L* is injected too deep or intramuscularly, the duration of the implant will be shorter because of a higher hyaluronic acid turnover rate. An injection too superficial may give blanching effects and bumps on the treatment site. If blanching is observed, i.e. the overlying skin turns a whitish colour, the injection should be stopped at once and the area massaged until it returns to a normal colour. Before injecting, press the rod carefully until a small droplet is visible at the tip of the needle.

Note! The extrusion force of the 29 G needle is approximately 50% less than with the 30 G needle for the 0.5 mL and 1 mL fill size. The reduced extrusion force should be considered when injecting the product.

The injection technique with regard to the depth of injection and the administered quantity may vary. The linear threading technique can be used to carefully lift up the wrinkle, but some physicians prefer a series of punctual injections or a combination of the two. During injection it is recommended that the eye of the needle should face upwards. The contour of the needle should be visible but not the colour of it. Inject *Restylane-L* while pulling the needle slowly backwards. Injection should stop just before the needle is pulled out from the skin to prevent material from leaking out from the injection site. In the treatment of lips, an enhanced vermilion border as well as fullness and pouting can be obtained. Please consult Medicis Aesthetics Canada Ltd. for details. Defects should be fully corrected, but not overcorrected, at each treatment session. If the skin of the patient is very loose, it is recommended that *Restylane-L* be injected on two separate occasions. The correction site should be massaged to conform to the contour of the surrounding tissues. For each treatment site a maximum dosage of 2 mL per treatment session is recommended. If the treated area is swollen directly after the injection, melting ice can be applied on the site for a short period. After the first treatment, additional implantations of *Restylane-L* may be necessary to achieve the desired level of correction. Periodic follow-up injections help sustain the desired degree of correction.

Note! The correct injection technique is crucial for the final result of the treatment. *Restylane-L* is only intended to be administered by authorized personnel in accordance with local legislation.

The syringe, the needle and any unused material must be discarded directly after the treatment session.

How supplied

Restylane-L is supplied in a disposable glass syringe with a Luer-Lok fitting. *Restylane-L* is co-packed with a sterilized needle(s) as indicated on the carton, either 29 G TW x 1/2" or 30 G x 1/2". A patient record label is a part of the syringe label. Remove it by pulling the flap marked with three small arrows (see picture D). This label is to be attached to patient records to ensure traceability of the product. The contents of the syringe are sterile. The number of units per package and the volume contained in each syringe is as stated on the outer package.

STERILE NEEDLES

- Follow national, local or institutional guidelines for use and disposal of medical sharp devices. Obtain prompt medical attention if injury occurs.
- To help avoid needle breakage, do not attempt to straighten a bent needle. Discard it and complete the procedure with a replacement needle.
- Do not resheat used needles. Recapping by hand is a hazardous practice and should be avoided.
- Discard unshielded needles in approved sharps collectors.

Shelf life and storage

As indicated on package. Store up to 25°C (77°F). Protect from sunlight and freezing. Refrigeration not required.

Canadian Patent No. 2,226,488

Manufactured by

Q-Med AB,
Seminariégatan 21, SE-752 28 Uppsala, Sweden
Phone +46(0)18 474 90 00, Fax +46(0)18 474 90 01
http://www.q-med.com, e-mail: info@q-med.com

Manufactured for

Medicis Aesthetics Canada Ltd./LTÉE
Toronto, Ontario
CANADA, M5V 1Y6
1-800-661-3376

IF THE PACKAGE IS DAMAGED, DO NOT USE

Symbols on packaging

- Refer to instructions for use
- For single use
- Sterile. The contents of the syringe have been sterilized by using moist heat.

Restylane, *Perlane* and *Restylane Fine Lines* are registered trademarks of HA North America Sales AB.

All other trademarks are the property of their respective owners.

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Restylane, *Perlane* et *Restylane Fine Lines* sont des marques appartenant à HA North America Sales AB.